

Performance Evaluation of Determine™ HIV-1/2 Ag/Ab Combo in plasma and whole blood from early HIV-1 infections

Silvina Masciotra

Special Studies and Training Activity Lead
HIV Diagnostics and Incidence Team
Atlanta, March 22, 2016

FDA-approved Ag/Ab HIV-1/2 screening assays

Assay	Manufacturer	specimen type	p24 detection	System
ADVIA Centaur HIV Ag/Ab Combo	SIEMENS	serum	with Ab	ADVIA Centaur or Centaur XP
ARCHITECT HIV Ag/Ab Combo	Abbott	serum/plasma	with Ab	ARCHITECT
GS HIV Combo Ag/Ab EIA	Bio-Rad	serum/plasma	with Ab	Evolis/automated microplate or manual
Bio-Plex 2200 HIV Ag-Ab	Bio-Rad	serum/plasma	differentiation	BioPlex 2200

- ❑ Instrumented or plate-based formats
- ❑ from ~30 min to 3 hours
- ❑ Not approved for blood or plasma donor screening

Alere Determine HIV-1/2 Ag/Ab Combo (DC)


- ❑ Manufacturer: Alere Scarborough, Inc.
- ❑ *In vitro*, visually read, rapid qualitative immunoassay
- ❑ Detection of HIV-1 p24 antigen and HIV-1/HIV-2 antibodies
- ❑ Point-of-care test to aid in the diagnosis of infection with HIV-1 and HIV-2, including an acute HIV-1 infection
- ❑ It may distinguish acute HIV-1 infection from established HIV-1 infection when the specimen is positive for HIV-1 p24 antigen and negative for anti-HIV-1 and anti-HIV-2 antibodies*

Alere Determine™ HIV-1/2 Ag/Ab Combo (DC)

- ❑ Human serum, plasma, and capillary (fingerstick) or venipuncture (venous) whole blood
- ❑ CLIA-waived for fingerstick whole blood
- ❑ 20 min using 50 µl of plasma/serum/whole blood



DC and CDC/APHL diagnostic algorithm



- ❑ Data on performance of DC as a screening test is limited
 - DC detected significantly fewer HIV-1 infections (52.6%) in early stages of seroconversion (38 vs. 20, $p < 0.0001$), but showed no difference with established infections
- ❑ Results from performance of DC on whole blood from individuals during early infections are inconsistent

Objectives

- ❑ To further evaluate the performance of DC as a screening assay in a subset of plasma specimens in the context of the new diagnostic algorithm
- ❑ To compare the performance of DC in simulated whole blood to plasma samples from HIV-1 seroconverters

STOP sample set

- 329 plasma specimens collected in San Francisco
 - Screening Targeted Populations to Interrupt On-going Chains of HIV Transmission with Enhanced Partner Notification (STOP)
 - Multi-site, prospective study evaluating methods to detect acute HIV-1 infections
 - Individuals were tested with FS whole blood with StatPak
 - Plasma specimens were tested with:
 - Abbott ARCHITECT, Bio-Rad Multispot and Abbott m2000 HIV-1 RNA viral load in San Francisco
 - DC at CDC and results compared to previous testing (McNemar's test)
 - Architect-false reactive plasma specimens were tested with Bio-Rad GS HIV Combo Ag/Ab EIA

HIV-1 seroconverters from the US

- 107 selected sequential plasma from 20 commercial seroconverters were used to simulate whole blood samples (40% hematocrit)
 - Whole blood samples were tested with DC and compared to results from plasma samples
 - DC-reactivity was calculated relative to days after the first Aptima-positive for both sample types

Results


DC performance in STOP specimens

Results from previous testing	total	HIV-1 RNA cop/ml median	Determine Combo results				DC positivity
			Ag-/Ab-	Ag+/Ab-	Ag+/Ab+	Ag-/Ab+	%
ARC-negative/HIV-1 RNA-positive	9	3.49E+03	0	0	0	0	0
ARC-positive/MS-negative/HIV-1 RNA-positive	74	1.65E+06	36	25	8	5	51.4
ARC-positive/MS-indeterminate/HIV-1 RNA-positive	11*	2.33E+06	3	1	3	4	72.7
ARC-positive/MS-positive or-undifferentiated	203**	1.47E+05	0	0	6	197	100
ARC-positive/HIV-1 RNA-negative	32	target not detected	0	0	0	0	0

* 10 and **30 VL results; ARC: Abbott ARCHITECT HIV Ag/Ab; MS: Bio-Rad Multispot HIV-1/HIV-2

- DC detected **54.1% (46/85)** of early HIV-1 infections
 - Difference was significantly different $p < 0.0001$
- All established HIV-1 infections were detected by DC
- No reactivity was observed among ARC-false reactive
 - All plasma were non-reactive with GS HIV Combo Ag/Ab EIA


Ag detection among 288 ARC-positive HIV-1 plasma specimens



* 15 of 17 Ag+Ab+ had VL results; $10^7 \Rightarrow >10^7$ cop/ml; 40 = <40 cop/ml; 0 = Target not detected

- Among 85 early HIV-1 infections, Ag detection was:
 - 33 of 75 (44.0%) of ARC-positive/Multispot-negative
 - Median VL 4.7×10^6 cop/ml
 - 4 of 11 (23.4%) of ARC-positive/Multispot-indeterminate
 - VL $>10 \times 10^7$ cop/ml in 3 samples


Marked reduction of Ag detection in early HIV-1 infections with whole blood compared to plasma



* One whole blood specimen was repeatedly invalid

- Overall DC reactivity in early HIV-1 infections was 91.1% with plasma and 56.4% with whole blood ($p < 0.0001$)

Nine of 20 SC showed delayed reactivity in whole blood compared to plasma



- As with plasma, DC with whole blood did not react in 8 samples collected before HIV-1 RNA was detected or in 14 RNA-positive only samples

DC performance with whole blood

- ❑ **Initial invalid results were obtained in 17.8% of whole blood specimens**
 - ❑ Two specimens that remained invalid after repeat testing were part of a seroconversion panel that never became positive up to 28 days after the first available Aptima positive (no further time points were tested)
- ❑ **Eight SC showed a median delay between plasma and whole blood of 6 days**
 - One SC also exhibited a second negative phase with whole blood for 7 days after being Ag+
- ❑ **Among 19 SC DC-plasma, the median time of DC reactivity since first Aptima positive**
 - Plasma: 6 days
 - Whole blood: 7 days

Conclusions

- ❑ DC used with plasma detected fewer specimens with early HIV-1 infection compared to an instrumented lab-based Ag/Ab assay
- ❑ DC with whole blood showed more invalid results than with plasma
- ❑ Results indicate that better sensitivity with plasma may be partially due to Ag detection

Considerations

- ❑ In this study, simulated whole blood was prepared from commercial seroconversion panels. Further studies will address sensitivity and specificity issues with FS whole blood
- ❑ In settings where lab-based testing is not feasible, DC might represent an advantage in detecting HIV-1 infections earlier
- ❑ More results are needed to evaluate the performance in the diagnostic algorithm to detect acute HIV-1 infections

Acknowledgments

CDC HIV lab

Wei Luo

Sarah Adams

Tara Smith

S. Michele Owen

STOP team

CDC

Philip J. Peters

New York City Department of Health and
Mental Hygiene

Emily Westheimer

San Francisco Department of Public Health

Stephanie Cohen

University of North Carolina

Cindy L. Gay

ICF International

Laura Hall

Thank you!

1600 Clifton Road NE, Atlanta, GA 30333

Telephone: 404-639-1004

svm6@cdc.gov

Visit: www.cdc.gov | Contact CDC at: 1-800-CDC-INFO or www.cdc.gov/info

The findings and conclusions in this report are those of the authors and do not necessarily represent the official position of the Centers for Disease Control and Prevention.

National Center for HIV/AIDS, Viral Hepatitis, STD, and TB Prevention
Division of HIV/AIDS Prevention

